

## Product Datasheet

### Mouse CD163 ELISA Kit (orb408595)

**Description** Mouse Scavenger receptor cysteine-rich type 1 protein M130(CD163) ELISA kit

**Reactivity** Mouse

**Range** 0.93 ng/ml-60ng/ml

**Target** CD163

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma, tissue homogenates

**Uniprot ID** [Q2VLH6](#)

**Sensitivity** 0.23 ng/ml

**Expiration Date** Please enquire.