

## Product Datasheet

### Mouse SFRP4 ELISA Kit (orb408601)

<b>Description</b>	Mouse Secreted frizzled-related protein 4(SFRP4) ELISA kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.47 ng/mL-30 ng/mL
<b>Target</b>	SFRP4
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, cell culture supernates, tissue homogenates
<b>Uniprot ID</b>	<a href="#">Q9Z1N6</a>
<b>Sensitivity</b>	0.11 ng/mL
<b>Expiration Date</b>	Please enquire.