

## Product Datasheet

### Mouse SORT1 ELISA Kit (orb408650)

**Description** Mouse Sortilin(SORT1) ELISA kit

**Reactivity** Mouse

**Range** 6.25 pg/mL-400 pg/mL

**Target** SORT1

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma, tissue homogenates

**Uniprot ID** [Q6PHU5](#)

**Sensitivity** 1.56 pg/mL

**Expiration Date** Please enquire.