

## Product Datasheet

### Mouse TNFSF13 ELISA Kit (orb408725)

<b>Description</b>	Mouse Tumor necrosis factor ligand superfamily member 13(TNFSF13) ELISA kit
<b>Reactivity</b>	Mouse
<b>Range</b>	93.75 pg/mL-6000 pg/mL
<b>Target</b>	TNFSF13
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, cell culture supernates, tissue homogenates, cell lysates
<b>Uniprot ID</b>	<a href="#">Q9D777</a>
<b>Sensitivity</b>	23.43 pg/mL
<b>Expiration Date</b>	Please enquire.