

Product Datasheet

Mouse TNFRSF17 ELISA Kit (orb408733)

Description	Mouse Tumor necrosis factor receptor superfamily member 17(TNFRSF17) ELISA kit
Reactivity	Mouse
Range	0.39 ng/ml- 20ng/ml
Target	TNFRSF17
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O88472
Sensitivity	0.009 ng/ml
Expiration Date	Please enquire.