

## Product Datasheet

### Mouse CXCL1 ELISA Kit (orb408819)

<b>Description</b>	Mouse chemokine (C-X-C motif) ligand 1 (melanoma growth stimulating activity, alpha) (CXCL1) ELISA k...
<b>Reactivity</b>	Mouse
<b>Range</b>	7.8 pg/mL-500 pg/mL
<b>Target</b>	CXCL1
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">P12850</a>
<b>Sensitivity</b>	1.95 pg/mL
<b>Expiration Date</b>	Please enquire.