

Product Datasheet

Mouse C-Peptide ELISA Kit (orb408824)

| | |
|--------------------------|------------------------------|
| Description | Mouse C-Peptide ELISA Kit |
| Reactivity | Mouse |
| Range | 2.86 ng/mL-50 ng/mL |
| Target | C-Peptide |
| Note | For research use only |
| Application notes | Detection Wavelength: 450 nm |
| Sample Types | serum, plasma |
| Sensitivity | 1.786 ng/mL |
| Expiration Date | Please enquire. |