

## Product Datasheet

### Mouse IL1RA ELISA Kit (orb408933)

|                          |   |
|--------------------------|---|
| <b>Description</b>       | Mouse Interleukin 1 receptor antagonist, IL-1ra ELISA kit |
| <b>Reactivity</b>        | Mouse   |
| <b>Range</b>             | 62.5 pg/mL-4000 pg/mL                                     |
| <b>Target</b>            | IL1RN   |
| <b>Note</b>              | For research use only                                     |
| <b>Application notes</b> | Detection Wavelength: 450 nm                              |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates                         |
| <b>Uniprot ID</b>        | <b>P25085</b>   |
| <b>Sensitivity</b>       | 15.6 pg/mL  |
| <b>Expiration Date</b>   | Please enquire.   |