

## Product Datasheet

# Mouse Complement fragment 5a ELISA Kit (orb408980)

<b>Description</b>	Mouse Complement fragment 5a,C5a ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	39 pg/mL-2500 pg/mL
<b>Target</b>	C5a
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Sensitivity</b>	9.75 pg/mL
<b>Expiration Date</b>	Please enquire.