

Product Datasheet

Mouse FGF15 ELISA Kit (orb408990)

Description	Mouse Fibroblast growth factor 15(FGF15) ELISA kit
Reactivity	Mouse
Range	15.6 pg/ml - 1000 pg/ml
Target	FGF15
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	O35622
Sensitivity	3.9 pg/ml
Expiration Date	Please enquire.