

Product Datasheet

Mouse TSH ELISA Kit (orb409000)

Description	Mouse Thyroid Stimulating Hormone,TSH ELISA Kit
Reactivity	Mouse
Range	0.3 µU/ml-12 µU/ml
Target	TSH
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates.
Uniprot ID	P12656
Sensitivity	0.15 µU/ml
Expiration Date	Please enquire.