

Product Datasheet

mouse ACTH ELISA Kit (orb409003)

| | |
|--------------------------|---|
| Description | Mouse adrenocorticotrophic hormone (ACTH) ELISA kit |
| Reactivity | Mouse |
| Range | 18.75 pg/mL-1200 pg/mL |
| Target | ACTH |
| Note | For research use only |
| Application notes | Detection Wavelength: 450 nm |
| Sample Types | serum, plasma |
| Sensitivity | 4.68 pg/mL |
| Expiration Date | Please enquire. |