

## Product Datasheet

### Mouse GNRH1 ELISA Kit (orb409006)

<b>Description</b>	Mouse gonadotropin-releasing hormone, GnRH ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	7.8 pg/ml - 500 pg/ml
<b>Target</b>	GNRH1
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">P13562</a>
<b>Sensitivity</b>	1.95 pg/ml
<b>Expiration Date</b>	Please enquire.