

Product Datasheet

Mouse CCL8 ELISA Kit (orb409011)

Description	Mouse monocyte chemotactic protein 2, MCP-2 ELISA kit
Reactivity	Mouse
Range	6.25 pg/mL-400 pg/mL
Target	CCL8
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9Z121
Sensitivity	1.56 pg/mL
Expiration Date	Please enquire.