

## Product Datasheet

### Mouse CCK ELISA Kit (orb409015)

<b>Description</b>	Mouse Cholecystokinin,CCK ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	31.25 pg/mL - 2000 pg/mL
<b>Target</b>	CCK
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>P09240</b>
<b>Sensitivity</b>	7.8 pg/mL
<b>Expiration Date</b>	Please enquire.