

Product Datasheet

Mouse SP-D ELISA Kit (orb409045)

Description	Mouse pulmonary surfactant-associated protein D (SP-D) ELISA kit
Reactivity	Mouse
Range	1.56 ng/mL-100 ng/mL
Target	SP-D
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P50404
Sensitivity	0.39 ng/mL
Expiration Date	Please enquire.