

## Product Datasheet

### Mouse SDF1 ELISA Kit (orb409142)

<b>Description</b>	Mouse Stromal cell derived factor 1 $\beta$ ,SDF-1 $\beta$ ELISA kit
<b>Reactivity</b>	Mouse
<b>Range</b>	31.25 pg/mL-2000 pg/mL
<b>Target</b>	CXCL12
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, cell culture supernates, tissue homogenates
<b>Uniprot ID</b>	<b>P40224</b>
<b>Sensitivity</b>	7.81 pg/mL
<b>Expiration Date</b>	Please enquire.