

Product Datasheet

Mouse FGF2 ELISA Kit (orb409154)

Description	Mouse basic fibroblast growth factor,bFGF ELISA Kit
Reactivity	Mouse
Range	3.12 pg/mL-200 pg/mL
Target	FGF2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	P15655
Sensitivity	0.78 pg/mL
Expiration Date	Please enquire.