

## Product Datasheet

### Mouse COL18A1/ES ELISA Kit (orb409159)

<b>Description</b>	Mouse Collagen alpha-1(XVIII) chain (COL18A1/ES) ELISA kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.156 ng/mL-10 ng/mL
<b>Target</b>	COL18A1/ES
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, cell culture supernates, tissue homogenates
<b>Uniprot ID</b>	<b>P39061</b>
<b>Sensitivity</b>	0.039 ng/mL
<b>Expiration Date</b>	Please enquire.