

Product Datasheet

Mouse CCL3 ELISA Kit (orb409175)

Description Mouse Macrophage Inflammatory Protein
1 α ,MIP-1 α ELISA kit

Reactivity Mouse

Range 7.8 pg/ml - 500 pg/ml

Target CCL3

Note For research use only

Application notes Detection Wavelength: 450 nm

Sample Types serum, plasma, tissue homogenates

Uniprot ID [P10855](#)

Sensitivity 1.95 pg/ml

Expiration Date Please enquire.