

Product Datasheet

Mouse MIP3 beta ELISA Kit (orb409179)

Description	Mouse macrophage inflammatory protein 3 β , MIP-3 β ELISA kit
Reactivity	Mouse
Range	39 pg/ml - 2500 pg/ml
Target	CCL19
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O70460
Sensitivity	9.75 pg/ml
Expiration Date	Please enquire.