

Product Datasheet

Mouse CD40LG ELISA Kit (orb409206)

Description	Mouse Soluble Cluster of differentiation 40 ligand,sCD40L ELISA Kit
Reactivity	Mouse
Range	0.274 ng/ml-200 ng/ml.
Target	CD40LG
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell culture supernates.
Uniprot ID	P27548
Sensitivity	0.274 ng/ml.
Expiration Date	Please enquire.