

Product Datasheet

Mouse TLR4 ELISA Kit (orb409210)

Description	Mouse Soluble Toll-like Receptor 4, sTLR4 ELISA Kit
Reactivity	Mouse
Range	0.156 ng/mL-10 ng/mL
Target	TLR4
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9QUK6
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.