

## Product Datasheet

### Mouse AMH ELISA Kit (orb409257)

**Description** Mouse anti-Mullerian hormone (AMH) ELISA kit

**Reactivity** Mouse

**Range** 0.4 ng/mL-14 ng/mL

**Target** AMH

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma

**Uniprot ID** [P27106](#)

**Sensitivity** 0.25 ng/mL

**Expiration Date** Please enquire.