

Product Datasheet

Mouse TCC C5b9 ELISA Kit (orb409260)

Description	Mouse terminal complement complex C5b-9,TCC C5b-9 ELISA Kit
Reactivity	Mouse
Range	0.78 ng/mL-50 ng/mL
Target	TCC C5b9
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Sensitivity	0.195 ng/mL
Expiration Date	Please enquire.