

## Product Datasheet

### Mouse VWF ELISA Kit (orb409286)

**Description** Mouse coagulation factor VIII(F VIII)  
ELISA Kit

**Reactivity** Mouse

**Range** 0.312 ng/mL-20 ng/mL

**Target** F8

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma

**Uniprot ID** [Q06194](#)

**Sensitivity** 0.078 ng/mL

**Expiration Date** Please enquire.