

## Product Datasheet

### Mouse PROK2 ELISA Kit (orb409302)

<b>Description</b>	Mouse Prokineticin-2(PROK2) ELISA kit
<b>Reactivity</b>	Mouse
<b>Range</b>	0.68 ng/ml - 500 ng/ml
<b>Target</b>	PROK2
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates.
<b>Uniprot ID</b>	<a href="#">Q9QXU7</a>
<b>Sensitivity</b>	0.68 ng/ml
<b>Expiration Date</b>	Please enquire.