

## Product Datasheet

### Mouse Lactoferrin ELISA Kit (orb409330)

<b>Description</b>	Mouse Lactoferrin,LTF/LF ELISA kit
<b>Reactivity</b>	Mouse
<b>Range</b>	Request Information
<b>Target</b>	LTF
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, breast milk
<b>Uniprot ID</b>	<b>P08071</b>
<b>Sensitivity</b>	Request Information
<b>Expiration Date</b>	Please enquire.