

Product Datasheet

Mouse CXCL5 ELISA Kit (orb409336)

Description	Mouse Epithelial neutrophil activating peptide 78,ENA-78 ELISA Kit
Reactivity	Mouse
Range	62.5 pg/ml - 4000 pg/ml
Target	CXCL5
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P50228
Sensitivity	15.6 pg/ml
Expiration Date	Please enquire.