

## Product Datasheet

### Mouse Adrenomedullin ELISA Kit (orb409352)

<b>Description</b>	Mouse Adrenomedullin ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	4.7 pg/mL-300 pg/mL
<b>Target</b>	ADM
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, cell culture supernates, tissue homogenates
<b>Uniprot ID</b>	<b>P97297</b>
<b>Sensitivity</b>	1.2 pg/mL
<b>Expiration Date</b>	Please enquire.