

## Product Datasheet

### Mouse CTLA4 ELISA Kit (orb409422)

**Description** Mouse cytotoxic T lymphocyte associated antigen 4,CTLA-4 ELISA Kit

**Reactivity** Mouse

**Range** 0.312 ng/mL-20 ng/mL

**Target** CTLA4

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma, tissue homogenates

**Uniprot ID** [P09793](#)

**Sensitivity** 0.078 ng/mL

**Expiration Date** Please enquire.