

Product Datasheet

Mouse ACE2 ELISA Kit (orb409453)

Description	Mouse Angiotensin converting enzyme 2(ACE2)ELISA Kit
Reactivity	Mouse
Range	6.25 pg/mL-400 pg/mL
Target	ACE2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q8R0I0
Sensitivity	1.56 pg/mL
Expiration Date	Please enquire.