

Product Datasheet

Mouse VEGFC ELISA Kit (orb409455)

Description	Mouse Vascular Endothelial cell Growth Factor C, VEGF-C ELISA kit
Reactivity	Mouse
Range	1.56 pg/mL-100 pg/mL
Target	VEGFC
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates, cell lysates
Uniprot ID	P97953
Sensitivity	0.39 pg/mL
Expiration Date	Please enquire.