

Product Datasheet

Mouse ANGPT2 ELISA Kit (orb409461)

Description	Mouse Angiopoietin 2,ANG-2 ELISA Kit
Reactivity	Mouse
Range	0.156 pg/mL-10 pg/mL
Target	ANGPT2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	O35608
Sensitivity	0.039 pg/mL
Expiration Date	Please enquire.