

Product Datasheet

Mouse ANGPTL4 ELISA Kit (orb409463)

| | |
|--------------------------|---|
| Description | Mouse Angiopoietin-related protein 4(ANGPTL4) ELISA kit |
| Reactivity | Mouse |
| Range | 3.12 pg/mL-200 pg/mL |
| Target | ANGPTL4 |
| Note | For research use only |
| Application notes | Detection Wavelength: 450 nm |
| Sample Types | serum, plasma, tissue homogenates |
| Uniprot ID | Q9Z1P8 |
| Sensitivity | 0.78 pg/mL |
| Expiration Date | Please enquire. |