

Product Datasheet

Mouse IGFBP1 ELISA Kit (orb409495)

Description Mouse Insulin-like growth factor binding protein 1,IGFBP-1 ELISA Kit

Reactivity Mouse

Range 0.312 ng/mL-20 ng/mL

Target IGFBP1

Note For research use only

Application notes Detection Wavelength: 450 nm

Sample Types serum, plasma, tissue homogenates

Uniprot ID [P47876](#)

Sensitivity 0.078 ng/mL

Expiration Date Please enquire.