

Product Datasheet

Mouse Neutrophil Elastase ELISA Kit (orb409544)

Description	Mouse neutrophil elastase,NE ELISA Kit
Reactivity	Mouse
Range	0.312 ng/mL-20 ng/mL
Target	ELA2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q3UP87
Sensitivity	0.078 ng/mL
Expiration Date	Please enquire.