

## Product Datasheet

### Pig IL1B ELISA Kit (orb409698)

<b>Description</b>	Pig Interleukin 1 $\beta$ ,IL-1 $\beta$ ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	3.9 pg/mL-250 pg/mL
<b>Target</b>	IL1B
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, cell culture supernates, tissue homogenates
<b>Uniprot ID</b>	<b>P26889</b>
<b>Sensitivity</b>	0.975 pg/mL
<b>Expiration Date</b>	Please enquire.