


Product Datasheet

MOCS1 antibody (orb40993)

Description	MOCS1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Molybdenum cofactor biosynthesis protein 1 protein (90-360AA)
Target	MOCS1
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NZB8
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat whole ce...