


Product Datasheet

MSRB3 antibody (orb40996)

Description	MSRB3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Methionine-R-sulfoxide reductase B3 protein (1-185AA)
Target	MSRB3
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8IXL7
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse large int...