

Product Datasheet

Rabbit Luteinizing hormone ELISA Kit (orb410016)

Description	Rabbit luteinizing hormone(LH)ELISA kit
Reactivity	Rabbit
Range	1.56 mIU/mL-100 mIU/mL
Target	LH
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q6IY74
Sensitivity	0.39 mIU/mL
Expiration Date	Please enquire.