

## Product Datasheet

### rabbit Immunoglobulin M ELISA Kit (orb410030)

<b>Description</b>	rabbit Immunoglobulin M,IgM ELISA Kit
<b>Reactivity</b>	Rabbit
<b>Range</b>	15.6 ng/mL-4000 ng/mL
<b>Target</b>	IgM
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma
<b>Sensitivity</b>	4 ng/mL
<b>Expiration Date</b>	Please enquire.