

Product Datasheet

Rabbit IGF1 ELISA Kit (orb410053)

Description	rabbit Insulin-like growth factor 1,IGF-1 ELISA Kit
Reactivity	Rabbit
Range	0.156 ng/mL-10 ng/mL
Target	IGF1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q95222
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.