

## Product Datasheet

### Rat CX3CR1 ELISA Kit (orb410208)

<b>Description</b>	Rat CX3C chemokine receptor 1(CX3CR1) ELISA kit
<b>Reactivity</b>	Rat
<b>Range</b>	23.5 pg/mL-1500 pg/mL
<b>Target</b>	CX3CR1
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates
<b>Uniprot ID</b>	<b>P35411</b>
<b>Sensitivity</b>	5.8 pg/mL
<b>Expiration Date</b>	Please enquire.