

## Product Datasheet

### Rat CYR61 ELISA Kit (orb410460)

<b>Description</b>	Rat Protein CYR61(CYR61) ELISA kit
<b>Reactivity</b>	Rat
<b>Range</b>	31.25 pg/mL-2000 pg/mL
<b>Target</b>	CYR61
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<a href="#">Q9ES72</a>
<b>Sensitivity</b>	7.81 pg/mL
<b>Expiration Date</b>	Please enquire.