

## Product Datasheet

### Rat CYP27B1 ELISA Kit (orb410557)

<b>Description</b>	Rat Vitamin D3 receptor(VDR) ELISA kit
<b>Reactivity</b>	Rat
<b>Range</b>	31.2 pg/mL-2000 pg/mL
<b>Target</b>	VDR
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>P13053</b>
<b>Sensitivity</b>	7.8 pg/mL
<b>Expiration Date</b>	Please enquire.