

## Product Datasheet

### Rat 8-OHdG ELISA Kit (orb410579)

<b>Description</b>	Rat 8-Hydroxy-desoxyguanosine,8-OHdG ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	0.312 ng/mL-20 ng/mL
<b>Target</b>	8-OHdG
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, cell culture supernates, tissue homogenates
<b>Sensitivity</b>	0.078 ng/mL
<b>Expiration Date</b>	Please enquire.