

## Product Datasheet

### Rat CTX-II ELISA Kit (orb410588)

<b>Description</b>	Rat C-telopeptide of type II collagen, CTX-II ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	0.312 ng/mL-20 ng/mL
<b>Target</b>	CTX-II
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Sensitivity</b>	0.078 ng/mL
<b>Expiration Date</b>	Please enquire.