

## Product Datasheet

### Rat TLR4 ELISA Kit (orb410671)

**Description** Rat toll-like receptor 4 (TLR4) ELISA Kit

**Reactivity** Rat

**Range** 0.625 ng/mL-40 ng/mL

**Target** TLR4

**Note** For research use only

**Application notes** Detection Wavelength: 450 nm

**Sample Types** serum, plasma, cell culture supernates,  
tissue homogenates

**Uniprot ID** [Q9QX05](#)

**Sensitivity** 0.156 ng/mL

**Expiration Date** Please enquire.