

## Product Datasheet

### Rat IFN gamma ELISA Kit (orb410698)

<b>Description</b>	Rat Interferon gamma; ,IFN-gamma; ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	0.625 pg/mL-40 pg/mL
<b>Target</b>	IFNG
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Wavelength: 450 nm
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>P01581</b>
<b>Sensitivity</b>	0.156 pg/mL
<b>Expiration Date</b>	Please enquire.